

The multi-site pDXA bone densitometer

Easy to measure multi-site scan with quick scanning speed

EXA-3000 is especially designed for scanning forearm and calcaneus of patient by DXA (Dual energy X-ray Absorptiometry) technology. It is fastest BMD (Bone Mineral Densitometry) in the world. The lifter, optional item it's specifically designed for EXA-3000 makes easy to lift up and down and it allows a patient to be more comfortable.

Features

- Quick scanning speed : 5sec
- Multi-site scan : Forearm, Calcaneus
- Precise measurement with DXA technology
- Easy operation
- Most comfortable design
- Low radiation
- Reliable result
- Trend report



Lifter (Optional)

Calcaneus

Forearm & Calcaneus report

The cutting edge design of pDXA bone densitometer

The art of beauty & gorgeous design of peripheral DXA system

EXA-PRESTO is portable pDXA (Peripheral Dual energy X-ray Absorptiometry) with unique and compact design and it performs quick measuring speed and more precise measurement with specially designed positioner with low radiation level and it gives you easy operation with conventional work station as well as laptop computer.

Features

- Quick scanning speed : 5sec.
- Scan site : Forearm
- Precise measurement with DXA technology
- Easy operation
- Most comfortable design
- Low radiation
- Reliable result
- Trend report
- 3D image analysis
- FRAX



Trend view

Forearm report

